


Shuttle-Dor



Features

+ Anteriorly attached tension band provides **dorsiflexion assist** during swing phase and active support and dampening at heel strike

Segment attaches to the Shuttle ankle joint stirrup



Standard


Available **separately or as combi-product** together with the Shuttle ankle joint

Precise assembly calibrated to 1/100mm for smooth operation free of play


Serial number for traceability

Clearly illustrated instructions for use

Fig. left



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product	AFO/KAFO	HKAFO/RGO	Article No.
Shuttle-Dor-22	up to 15 kg	up to 15 kg	SH-Dor-22
Shuttle-Dor-25	up to 25 kg	up to 25 kg	SH-Dor-25
Shuttle-Dor-28	up to 35 kg	up to 35 kg	SH-Dor-28
Shuttle-Dor-30	up to 45 kg	up to 45 kg	SH-Dor-30
Shuttle-Dor-32	up to 55 kg	up to 55 kg	SH-Dor-32
Shuttle-Dor-34	(up to 70 kg) not yet available		SH-Dor-34

