

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics


Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026

Addition Separator-E

Simple Operation in 3 Steps






Self-contained quick-release system for separating **KAFO** into **K+AFO**, **HKAFO** into **HK+AFO**, or **RGO** into **AFO**.
Can be used with all orthotic knee and ankle joints.



Weight Classification (body weight)

Measurements

Article No.	AFO/KAFO	Width	Thickness	Length
Separator-E-klein-16	up to 30 kg	16 mm	12,0 mm	80 mm
Separator-E-mittel-19	up to 70 kg	19 mm	11,5 mm	90 mm
Separator-E-groß-22	up to 110 kg	22 mm	12,0 mm	110 mm



Modular

Detachable

Locked

No Grinding Necessary

Growable

Water Resistant



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com