

ORTHO
SYSTEMS

2000
2026
Years in motion



Excellence in Neuro-Orthopedics


Your Patient – Our Contribution

Catalogue
Excerpt

Orthotics 14.0
Valid from July 2026


On the Go with the Separator

Our Separators enable quick and easy separation, for example, of a KAFO into an AFO. This allows flexible use of the AFO in various everyday situations—such as in a car, on a train, in a wheelchair, or when cycling. In such cases, an orthosis extending above the knee can be restrictive and impractical. The orthosis can then be quickly reconnected at any time.



Donning, Doffing and Flexibility with the Separator

Donning an orthosis should be as simple as possible, as it is an everyday task. The more straightforward the process, the more naturally it can be integrated into everyday life. This not only increases independence, but also improves the quality of life for both users and caregivers. The ability to separate and reconnect an orthosis opens up a wide range of options in everyday life for the wearer. For example, the thigh section can be attached for **longer walking distances** yet selectively removed for **therapeutic purposes**. At the same time, the two-part design makes donning and doffing intuitive and user-friendly, whether performed independently or with **assistance**, while still allowing the orthosis to be reconnected quickly and securely in both cases.



KAFO for more knee stability

KAFO with Extension Assist –
usable as an AFO, e.g. for therapy

Assists donning while seated



Contact

You can reach us from
by telephone
per Fax
per E-Mail
Address
Contact persons

9 am to 5 pm (Germany)
+49 89 904 7575 0
+49 89 904 7575 5
kontakt@ortho-systems-mail.de
Merowingerstr. 14, 85551 Kirchheim, Germany
Jeff Matwichuk and Thomas Böckh

www.ortho-systems-inmotion.com